

GL GONIOSPECTROMETER MEASUREMENT TAG LIST

Report Wizard Sample Description Lines

Client name: *Azzardo Lighting*
Model name: *Lauren A 3-Step Dimm*
Item No: *AZ5258*
Serial No.:
Made by: *A. Nawrot*

Eulumdat Fields

Dimensions of luminaire

Company: *Azzardo*
Report number: *AZ5258*
LuminaireName: *Lauren A 3-Step Dimm GO*
LuminaireNumber: *AZ5258*
FileName: *AZ5258.Idt*
Date: *2024-06-19*

Dimensions of luminaire

Length: *1200 mm*
Width: *60 mm*
Height: *1500 mm*

Dimensions of luminous area

Length: *1200 mm*
Width: *60 mm*
C0-plane *0 mm*
C90-plane *0 mm*
C180-plane *0 mm*
C270-plane *0 mm*

Lamps:

Number of lamps: *-1*

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03

Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432

EUR: IBAN 06 1750 0012 0000 0000 2709 9467

SWIFT BIC: RCBWPLPW

NIP: 783-171-97-20

REGON 360317241

KRS: 0000537079

Kapitał 200.000PLN

Type of lamps: *LED*
Total luminous flux: *2073,75 lm*
Color temperatur: *2845 K*
Color rendering index: *81,14*
Wattage including ballast: *40,71 W*

Utilization factors:

k = 0.60	<i>0,235</i>
k = 0.80	<i>0,309</i>
k = 1.00	<i>0,375</i>
k = 1.25	<i>0,442</i>
k = 1.50	<i>0,496</i>
k = 2.00	<i>0,580</i>
k = 2.50	<i>0,632</i>
k = 3.00	<i>0,673</i>
k = 4.00	<i>0,730</i>
k = 5.00	<i>0,769</i>

Electrical parameters

<i>U</i>	<i>= 2,458 V</i>
<i>I</i>	<i>= 0,000 A</i>
<i>P</i>	<i>= 40,710 W</i>
<i>Power factor</i>	<i>= 0,000</i>

<i>Power apparent</i>	<i>= 0,000 VA</i>
<i>Power Reactive</i>	<i>= 0,000 var</i>
<i>Frequency Voltage</i>	<i>= 243,413 Hz</i>
<i>THD</i>	<i>= 0,000 %</i>

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432

EUR: IBAN 06 1750 0012 0000 0000 2709 9467

SWIFT BIC: RCBWPLPWNIP: 783-171-97-20

REGON 360317241

KRS: 0000537079

Kapitał 200.000PLN

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432

EUR: IBAN 06 1750 0012 0000 0000 2709 9467

SWIFT BIC: RCBWPLPWNIP: 783-171-97-20

REGON 360317241

KRS: 0000537079

Kapitał 200.000PLN

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432

EUR: IBAN 06 1750 0012 0000 0000 2709 9467

SWIFT BIC: RCBWPLPWNIP: 783-171-97-20

REGON 360317241

KRS: 0000537079


Kapitał 200.000PLN

Temperature = 25 C

Goniometric measurement results

Total flux (utilized luminarie flux)	2073,75 lm
Flux in lower hemisphere	100,00 %
Flux in upper hemisphere	-0,00 %
Maximum luminous intensity	556,55 cd
Light output ratio (LOR)	100,00 %
Luminous efficacy	50,94 lm/W

Light intensity distribution diagram



CIE 1931 diagram

x	0,4473
y	0,4056
Angular Colour Uniformity	0,0019

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432


EUR: IBAN 06 1750 0012 0000 0000 2709 9467

SWIFT BIC: RCBWPLPWNIP: 783-171-97-20

REGON 360317241

KRS: 0000537079

Kapitał 200.000PLN



Beam cone

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432


EUR: IBAN 06 1750 0012 0000 0000 2709 9467

SWIFT BIC: RCBWPLPWNIP: 783-171-97-20

REGON 360317241

KRS: 0000537079

Kapitał 200.000PLN



Stabilization

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432


EUR: IBAN 06 1750 0012 0000 0000 2709 9467

SWIFT BIC: RCBWPLPWNIP: 783-171-97-20

REGON 360317241

KRS: 0000537079

Kapitał 200.000PLN



3D Chart

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432


EUR: IBAN 06 1750 0012 0000 0000 2709 9467

SWIFT BIC: RCBWPLPWNIP: 783-171-97-20

REGON 360317241

KRS: 0000537079

Kapitał 200.000PLN



Legend

— C0-C180

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432

EUR: IBAN 06 1750 0012 0000 0000 2709 9467

SWIFT BIC: RCBWPLPWNIP: 783-171-97-20

REGON 360317241

KRS: 0000537079

Kapitał 200.000PLN



- C90-C270
- G0-G180

U(%f)

I(%f)

**Ext
end
ed
elec
tric
al
rep
ort**

P(%f)			
-------	--	--	--

0	195,970	228,252	-481,919
1	100,000	100,000	100,000
2	38,147	0	0
3	11,475	0	0
4	13,115	0	0
5	9,725	0	0
6	9,404	0	0
7	6,879	0	0
8	5,243	0	0
9	6,134	0	0
10	4,299	0	0
11	3,854	0	0
12	3,652	0	0
13	3,075	0	0
14	2,856	0	0
15	2,831	0	0
16	2,289	0	0
17	2,121	0	0
18	2,379	0	0
19	1,879	0	0
20	1,826	0	0
21	2,005	0	0

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432

EUR: IBAN 06 1750 0012 0000 0000 2709 9467

SWIFT BIC: RCBWPLPWNIP: 783-171-97-20

REGON 360317241

KRS: 0000537079

Kapitał 200.000PLN

22	1,736	0	0
23	1,662	0	0
24	1,948	0	0
25	1,563	0	0
26	1,476	0	0
27	1,851	0	0
28	1,426	0	0
29	1,443	0	0
30	1,621	0	0
31	1,212	0	0
32	1,335	0	0
33	1,398	0	0
34	1,169	0	0
35	1,137	0	0
36	1,264	0	0
37	1,047	0	0
38	1,009	0	0
39	1,129	0	0
40	0,924	0	0
41	0,968	0	0
42	1,082	0	0
43	0,883	0	0
44	0,919	0	0
45	1,050	0	0
46	0,876	0	0
47	0,752	0	0
48	0,987	0	0
49	0,798	0	0
50	0,756	0	0

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432

EUR: IBAN 06 1750 0012 0000 0000 2709 9467

SWIFT BIC: RCBWPLPWNIP: 783-171-97-20

REGON 360317241

KRS: 0000537079

Kapitał 200.000PLN