

GL GONIOSPECTROMETER MEASUREMENT TAG LIST

Report Wizard Sample Description Lines

Client name: *Azzardo Lighting*
Model name: *AZ5365 + AZ 5366 - 3000K*
Item No: *AZ5365 + AZ 5366 - 3000K*
Serial No.:
Made by: *A. Nawrot*

Eulumdat Fields

Dimensions of luminaire

Company: *Azzardo*
Report number: *AZ5365 , AZ5366 - 3000*
LuminaireName: *Panca 30 CCT Switch DIM - 3000K*
LuminaireNumber: *AZ5365-AZ5366-3000K*
FileName: *AZ5365_AZ5366_3000.Idt*
Date: *2023-02-09*

Dimensions of luminaire

Length: *300 mm*
Width: *0 mm*
Height: *25 mm*

Dimensions of luminous area

Length: *280 mm*
Width: *0 mm*
C0-plane *0 mm*
C90-plane *0 mm*
C180-plane *0 mm*
C270-plane *0 mm*

Lamps:

Number of lamps: -1
Type of lamps:
Total luminous flux: 1917,93 lm
Color temperatur: 2934 K
Color rendering index: 89,31
Wattage including ballast: 24,28 W

Utilization factors:

k = 0.60	0,292
k = 0.80	0,377
k = 1.00	0,452
k = 1.25	0,525
k = 1.50	0,582
k = 2.00	0,670
k = 2.50	0,720
k = 3.00	0,758
k = 4.00	0,811
k = 5.00	0,844

Electrical parameters

<i>U</i>	= 229,632 V
<i>I</i>	= 0,111 A
<i>P</i>	= 24,280 W
<i>Power factor</i>	= 0,951


<i>Power apparent</i>	= 25,530 VA
<i>Power Reactive</i>	= -5,905 var
<i>Frequency Voltage</i>	= 50,002 Hz
<i>THD</i>	= 0,000 %

Temperature = 25 C

Goniometric measurement results


Total flux (utilized luminarie flux)	1917,93 <i>lm</i>
Flux in lower hemisphere	100,00 %
Flux in upper hemisphere	0,00 %
Maximum luminous intensity	653,06 <i>cd</i>
Light output ratio (LOR)	100,00 %
Luminous efficacy	78,99 <i>lm/W</i>

Light intensity distribution diagram



CIE 1931 diagram

x 0,4386
y 0,3994
Angular Colour Uniformity 0,0009



GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczkowo; Poland

Tel.: +48 61 819 40 03

Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432

EUR: IBAN 06 1750 0012 0000 0000 2709 9467

SWIFT BIC: RCBWPLPW

NIP: 783-171-97-20

REGON 360317241

KRS: 0000537079

Kapitał 200.000 PLN

Beam cone

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03

Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432

EUR: IBAN 06 1750 0012 0000 0000 2709 9467


SWIFT BIC: RCBWPLPW

NIP: 783-171-97-20


REGON 360317241

KRS: 0000537079

Kapitał 200.000PLN



Stabilization



3D Chart

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03

Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432

EUR: IBAN 06 1750 0012 0000 0000 2709 9467


SWIFT BIC: RCBWPLPW

NIP: 783-171-97-20

REGON 360317241

KRS: 0000537079

Kapitał 200.000 PLN



Legend

- - C0-C180
- - C90-C270
- - G0-G180

GL Optic R&D Center Sp. z o.o.

ul. Poznańska 70

PL-62-040 Puszczykowo; Poland

Tel.: +48 61 819 40 03

Raiffeisen Bank Polska S.A.

PLN: IBAN 78 1750 0012 0000 0000 2709 9432

EUR: IBAN 06 1750 0012 0000 0000 2709 9467

SWIFT BIC: RCBWPLPW

NIP: 783-171-97-20

REGON 360317241

KRS: 0000537079

Kapitał 200.000PLN

Extended electrical report

	U(%f)	I(%f)	P(%f)
0	0	0	0
1	0	0	0
2	0	0	0
3	0	0	0
4	0	0	0
5	0	0	0
6	0	0	0
7	0	0	0
8	0	0	0
9	0	0	0
10	0	0	0
11	0	0	0
12	0	0	0
13	0	0	0
14	0	0	0
15	0	0	0
16	0	0	0
17	0	0	0
18	0	0	0
19	0	0	0
20	0	0	0
21	0	0	0
22	0	0	0
23	0	0	0
24	0	0	0
25	0	0	0
26	0	0	0
27	0	0	0
28	0	0	0
29	0	0	0
30	0	0	0
31	0	0	0
32	0	0	0
33	0	0	0
34	0	0	0
35	0	0	0
36	0	0	0
37	0	0	0
38	0	0	0
39	0	0	0
40	0	0	0
41	0	0	0
42	0	0	0

GL Optic R&D Center Sp. z o.o.

Raiffeisen Bank Polska S.A.

NIP: 783-171-97-20

ul. Poznańska 70

PLN: IBAN 78 1750 0012 0000 0000 2709 9432

REGON 360317241

PL-62-040 Puszczykowo; Poland

EUR: IBAN 06 1750 0012 0000 0000 2709 9467

KRS: 0000537079

Tel.: +48 61 819 40 03

SWIFT BIC: RCBWPLPW

Kapitał 200.000PLN

43	0	0	0
44	0	0	0
45	0	0	0
46	0	0	0
47	0	0	0
48	0	0	0
49	0	0	0
50	0	0	0